

# IMAGING PROTOCOL – CT DCN CTA Head & Neck Siemen's Edge+128

- CTDI: 0-4yr: ≤ 15mGy 5-16: ≤25mGy
- WO and post scans only with radiologist approval.
- Setup: Supine, LAT scout from vertex through the level of the aortic arch, no gantry angle.

## PT Positioning:

- Supine, AP and Lateral scout, no gantry angle
- Scout should extend through the aortic arch for smart prep/bolus tracking
- Retract shoulders as much as possible & tilt the patient's head so that a line connecting the lateral canthus of the eye and the EAC is perpendicular to the CT tabletop.

DFOV: Preferred 15cm

## Scan Parameters:

- Scan from vertex to below the aortic arch (to include most of arch)
- Contrast:
  - 2cc per kg of 350 mg (OMNI 300 if >50 lbs{22.6kg}), not to exceed 100 ml unless ok by radiologist
  - Bolus Tracking in ascending aorta, start scanning upon entry of contrast at the level of the aortic arch or trigger at 40 HU.
  - 10 second monitoring delay on bolus tracking/monitoring

## PACS SERIES

1. SCOUT & Patient Protocol (Dose Report)
2. ST AXL 1 x 1
3. BONE AXL 1 x 1
4. MIP COR 1 X 1
5. MIP SAG 1 X 1
6. MIP LAO 1 x 1
7. MIP RAO 1 x 1

## Acquisition Parameters

0-40 lbs		40-110lb		110lb +/-Adult	
Scan Type	Spiral	Scan Type	Spiral	Scan Type	Spiral
Pitch	1.4	Pitch	1.4	Pitch	1.2
Detector Configuration	128 x 0.6	Detector Configuration	128 x 0.6	Detector Configuration	128 x 0.6
Slice Thickness	1	Slice Thickness	1	Slice Thickness	1
Rotation Time	0.33	Rotation Time	0.5	Rotation Time	0.28
Care Dose	on	Care Dose	on	Care Dose	on
Care kV	on	Care kV	on	Care kV	on
Dose Optimization	9	Dose Optimization	9	Dose Optimization	9
Quality Ref mAs	156	Quality Ref mAs	156	Quality Ref mAs	84
Ref kV	80	Ref kV	80	Ref kV	120

## Reconstruction Parameters

<b>Recon 1 ST Axial</b>	
<b>Kernel</b>	<b>Br38</b>
<b>ADMIRE</b>	<b>3</b>
<b>Window</b>	<b>Soft Tissue 400 40</b>
<b>Slice Thickness</b>	<b>1 x 1</b>
<b>Recon 2 Bone Axial</b>	
<b>Kernel</b>	<b>Br59</b>
<b>ADMIRE</b>	<b>3</b>
<b>Window</b>	<b>Shoulder 2000 x 400</b>
<b>Slice Thickness</b>	<b>1 x 1</b>
<b>Recon 3 MIP Coronal</b>	<b>Recon 4 MIP Sagittal</b>
<b>Kernel: Bv38</b>	<b>Kernel: Bv38</b>
<b>ADMIRE: 3</b>	<b>ADMIRE: 3</b>
<b>Window: CT Angio 700 80</b>	<b>Window: CT Angio 700 80</b>
<b>Slice Thickness 1 x 1</b>	<b>Slice Thickness 1 x 1</b>
<b>LAO and RAO 45 ° to see carotid split</b>	
<b>Recon 5 LAO MIP</b>	<b>Recon 6 RAO MIP</b>
<b>Kernel: Bv38</b>	<b>Kernel: Bv38</b>
<b>ADMIRE: 3</b>	<b>ADMIRE: 3</b>
<b>Window: CT Angio 700 80</b>	<b>Window: CT Angio 700 80</b>
<b>Slice Thickness: 1 x 1</b>	<b>Slice Thickness: 1 x 1</b>