

Ascension Seton GE Revolution 256 Trauma Body Post Contrast 09/01/2021

Labs should be obtained from the emergency department

Scan Parameters:

Intravenous contrast at the discretion of the Radiologist

Patient Positioning:

- Both arms raised above the head for optimal image quality.
- If only one arm can be raised, secure the opposing arm along the patient's side
- If both arms are unable to be raised, secure the arms along the patient's side

Post contrast- To include above apices through symphysis pubis

Delay- (*Delay time is dependent on the urgency induced by the patient condition*)
To include above diaphragm through bladder

CTDI: ≤ 25 mGy

Detector Rows: 128

Detector Configuration: 0.625

Rotation time: 0.35

Pitch: 0.992

mA: 550

kV: 120

Auto mA: on

Smart mA: on

Noise Index: 15.4

Approved by Ascension Seton Medical Physicist: 2021.09.02

Approved by ARA Radiologists: 2021.09.27

Reconstruction Parameters:

Post Contrast

Recon 1: Axial Soft Tissue through the entire scan range

2.5 mm slice thickness
2.5 mm slice increment
Algorithm: Standard
W/L: 400/40
ASIR: 50

Recon 2: Axial Lung through the entire scan range

2.5 mm slice thickness
2.5 mm slice increment
Algorithm: Lung
W/L: 1600/-600
ASIR: none

Recon 3: (reformat set)

Axial Bone Recon of the Chest used for 3D Ribs:

0.625 mm slice thickness
0.3125 mm slice increment
Algorithm: Bone
W/L: 400/40

* 3D Ribs – use recon 3 to create a 3D Volume Rendered rotation and 3D Bone Window rotation of the ribs. Each series will have approx. 24 images in a Left to Right rotation.

Recon 4: Coronal Body

3.0 mm slice thickness
3.0 mm slice increment
Algorithm: General
W/L: 400/40
ASIR: None

Include the following reconstructions when a Thoracic and Lumbar Spine are ordered

Axial Soft Tissue Thoracic Spine

1.25 mm slice thickness
1.25 mm slice increment
DFOV 20 cm
Algorithm: Standard
W/L 400/40
ASIR: 50%

Axial Bone Thoracic Spine

1.25 mm slice thickness
1.25 mm slice increment
DFOV 20 cm
Algorithm: Bone Plus
W/L: 3000/500
ASIR: 50%

Thoracic Spine Sagittal Coronal MPR

2.0 mm slice thickness
1.5 mm slice increment
Algorithm: Bone
W/L: 3000/500

Thoracic Spine Bone MPR

2.0 mm slice thickness
1.5 mm slice increment
Algorithm: Bone
W/L: 3000/500

Axial Soft Tissue Lumbar Spine

1.25 mm slice thickness
1.25 mm slice increment
DFOV 15 cm
Algorithm: Standard
W/L: 400/40
ASIR: 50%

Axial Bone Lumbar Spine

1.25 mm slice thickness
1.25 mm slice increment
DFOV 15 cm
Algorithm: Bone Plus
W/L: 3000/500
ASIR: 50%

Lumbar Spine Sagittal MPR

2.0 mm slice thickness
1.5 mm slice increment
Algorithm: Bone
W/L: 3000/500

Lumbar Spine Coronal MPR

2.0 mm slice thickness
1.5 mm slice increment
Algorithm: Bone
W/L: 3000/500

GE Trauma Body Post Contrast 08/26/2021

Labs should be obtained from the emergency department

Scan Parameters:

Intravenous contrast at the discretion of the Radiologist

Patient Positioning:

- Both arms raised above the head for optimal image quality.
- If only one arm can be raised, secure the opposing arm along the patient's side
- If both arms are unable to be raised, secure the arms along the patient's side

Post contrast- To include above apices through symphysis pubis

Delay- (*Delay time is dependent on the urgency induced by the patient condition*)

To include above diaphragm through bladder

(Below Parameters are defined by a GE Optima 660, please adjust your GE scanner accordingly)

CTDI: ≤ 25 mGy

Detector Rows: 64

Detector Configuration: 0.625

Rotation time: 0.5

Pitch: 0.984:1

mA: 440

kV: 120

Auto mA: on

Smart mA: on

Noise Index: 15

Approved by Ascension Seton Medical Physicist: 2021.09.02

Approved by ARA Radiologists: 2021.09.27

Reconstruction Parameters:

(Below Parameters are defined by a GE Optima 660, please adjust your scanner accordingly)

Post Contrast

Recon 1: Axial Soft Tissue through the entire scan range

5 mm slice thickness
5 mm slice increment
Algorithm: Standard
Window: Mediastinum
ASIR: 20

Recon 2: Axial Lung through the entire scan range

5 mm slice thickness
5 mm slice increment
Algorithm: Lung
Window: Lung
ASIR: none

Recon 3: (reformat set)

Axial Bone Recon of the Chest used for 3D Ribs:

1.25 mm slice thickness
0.625 mm slice increment
Algorithm: Bone
Window: Bone

* 3D Ribs – use recon 3 to create a 3D Volume Rendered rotation and 3D Bone Window rotation of the ribs. Each series will have approx. 24 images in a Left to Right rotation.

Recon 4: Coronal Body

2.5 mm slice ASIR: none
thickness
2.5 mm slice increment
Algorithm: Standard
Window: Mediastinum
ASIR: 20

Recon 5: Axial Delay

5 mm slice thickness
5 mm slice increment
Algorithm: Standard
Window: Mediastinum
ASIR: 20

Sternum Recon*

Sagittal MPR of the Chest

2.5 mm slice thickness
2.5 mm slice increment
Algorithm: Bone
Window: Bone
ASIR: none

Axial Spine Recon*

Without Order

Sagittal MPR entire spine
2.5 mm slice thickness
2.5 mm slice increment
Algorithm: Bone
Window: Bone
ASIR: none

Include the following reconstructions when a Thoracic and Lumbar Spine are ordered

Focused DFOV

Axial Soft Tissue Thoracic Spine

2.5 mm slice thickness
2.5 mm slice increment
DFOV 15 cm
Algorithm: Standard
Window: Mediastinum
ASIR: none

Axial Bone Thoracic Spine

2.5 mm slice thickness
2.5 mm slice increment
DFOV 15 cm
Algorithm: Bone
Window: Bone
ASIR: none

Thoracic Spine Sagittal Coronal MPR

1.25 mm slice thickness
1.25 mm slice increment
Algorithm: Bone
Window: Bone

Thoracic Spine Coronal MPR

1.25 mm slice thickness
1.25 mm slice increment
Algorithm: Bone
Window: Bone

Axial Soft Tissue Lumbar Spine

2.5 mm slice thickness
2.5 mm slice increment
DFOV 15 cm
Algorithm: Standard
Window: Mediastinum
ASIR: none

Axial Bone Lumbar Spine

2.5 mm slice thickness
2.5 mm slice increment
DFOV 15 cm
Algorithm: Bone
Window: Bone
ASIR: none

Lumbar Spine Sagittal MPR

1.25 mm slice thickness
1.25 mm slice increment
Algorithm: Bone
Window: Bone

Lumbar Spine Coronal MPR

1.25 mm slice thickness
1.25 mm slice increment
Algorithm: Bone
Window: Bone

Ascension Seton Siemens 128 Trauma Body Post Contrast

09/01/2021

Labs should be obtained from the emergency department

Scan Parameters:

Intravenous contrast at the discretion of the Radiologist

Patient Positioning:

- Both arms raised above the head for optimal image quality.
- If only one arm can be raised, secure the opposing arm along the patient's side
- If both arms are unable to be raised, secure the arms along the patient's side

Post contrast- To include above apices through symphysis pubis

Delay- (*Delay time is dependent on the urgency induced by the patient condition*)

To include above diaphragm through bladder

CTDI: \leq 25 mGy

Quality Reference mAs: 210

Pitch: 0.6

CARE Dose4D: ON

Detector Rows: 128

Detector Configuration 0.6

CARE kV: 120

Rotation time: 0.5

Approved by Ascension Seton Medical Physicist: 2021.09.02

Approved by ARA Radiologists: 2021.09.27

Reconstruction Parameters:

Post Contrast

Recon 1: Axial Soft Tissue through the entire scan range

5 mm slice thickness

5 mm slice increment

Kernel: Br38

Window: Mediastinum

SAFIRE: 2

Recon 2: Axial Lung through the entire scan range

5 mm slice thickness

5 mm slice increment

Kernel: Br59

Window: Lung

SAFIRE: 0

Recon 3: (reformat set)

Axial Bone Recon of the Chest used for 3D Ribs:

1.0 mm slice thickness

0.5 mm slice increment

Kernel: B70s Sharp

Window: Osteo

SAFIRE: 0

* 3D Ribs – use recon 3 to create a 3D Volume Rendered rotation and 3D Bone Window rotation of the ribs. Each series will have approx. 24 images in a Left to Right rotation.

Recon 4: Coronal Body

5 mm slice thickness

5 mm slice increment

Kernel: I41f medium

Window: Abdomen

SAFIRE: 2

Recon 5: Axial Delay

5 mm slice thickness

5 mm slice increment

Kernel: I41f medium

Window: Abdomen

SAFIRE: 2

Sternum Recon*

Sagittal MPR of the Chest

2 mm slice thickness

2 mm slice increment

Kernel: B70s Sharp

Window: Osteo

SAFIRE: 0

Axial Spine Recon*

Without Order

Sagittal MPR entire spine

2 mm slice thickness

2 mm slice increment

Kernel: B70s Sharp

Window: Osteo

SAFIRE: 0

Include the following reconstructions when a Thoracic and Lumbar Spine are ordered

Focused DFOV

Axial Soft Tissue Thoracic Spine

2.5 mm slice thickness

2.5 mm slice increment

DFOV 15 cm

Kernel: I41f medium

Window: Spine

SAFIRE: 2

Axial Bone Thoracic Spine

2.5 mm slice thickness

2.5 mm slice increment

DFOV 15 cm

Kernel: B70s Sharp

Window: Osteo

SAFIRE: 0

Thoracic Spine Sagittal Coronal MPR

1.25 mm slice thickness

1.25 mm slice increment

Kernel: B70s Sharp

Window: Osteo

SAFIRE: 0

Thoracic Spine Coronal MPR

1.25 mm slice thickness

1.25 mm slice increment

Kernel: B70s Sharp

Window: Osteo

SAFIRE: 0

Axial Soft Tissue Lumbar Spine

2.5 mm slice thickness

2.5 mm slice increment

DFOV 15 cm

Kernel: I41f medium

Window: Spine

SAFIRE: 2

Axial Bone Lumbar Spine

2.5 mm slice thickness

2.5 mm slice increment

DFOV 15 cm

Kernel: B70s Sharp

Window: Osteo

SAFIRE: 0

Lumbar Spine Sagittal MPR

1.25 mm slice thickness

1.25 mm slice increment

Kernel: B70s Sharp

Window: Osteo

SAFIRE: 0

Siemens Trauma Body Post Contrast 08/26/2021

Labs should be obtained from the emergency department

Scan Parameters:

Intravenous contrast at the discretion of the Radiologist

Patient Positioning:

- Both arms raised above the head for optimal image quality.
- If only one arm can be raised, secure the opposing arm along the patient's side
- If both arms are unable to be raised, secure the arms along the patient's side

Post contrast- To include above apices through symphysis pubis

Delay- (*Delay time is dependent on the urgency induced by the patient condition*)
To include above diaphragm through bladder

(Below Parameters are defined by a Siemens Definition AS 64, please adjust your Siemens scanner accordingly)

CTDI: ≤ 25 mGy

Quality Reference mAs: 200

Pitch: 1.0

CARE Dose4D: ON

Detector Rows: 64

Detector Configuration 0.6

CARE kV: 120 (Optimize slider 7)

Rotation time: 0.5

Approved by Ascension Seton Medical Physicist: 2021.09.02

Approved by ARA Radiologists: 2021.09.27

Reconstruction Parameters:

(Below Parameters are defined by a Definition AS 64, please adjust your scanner accordingly)

Post Contrast

Recon 1: Axial Soft Tissue through the entire scan range

5 mm slice thickness

5 mm slice increment

Kernel: I41f medium

Window: Abdomen

SAFIRE: 2

Recon 2: Axial Lung through the entire scan range

5 mm slice thickness

5 mm slice increment

Kernel: I71f very sharp ASA

Window: Lung

SAFIRE: 0

Recon 3: (reformat set)

Axial Bone Recon of the Chest used for 3D Ribs:

1.0 mm slice thickness

0.5 mm slice increment

Kernel: B70s Sharp

Window: Osteo

SAFIRE: 0

* 3D Ribs – use recon 3 to create a 3D Volume Rendered rotation and 3D Bone Window rotation of the ribs. Each series will have approx. 24 images in a Left to Right rotation.

Recon 4: Coronal Body

2 mm slice thickness

2 mm slice increment

Kernel: I41f medium

Window: Abdomen

SAFIRE: 2

Recon 5: Axial Delay

5 mm slice thickness

5 mm slice increment

Kernel: I41f medium

Window: Abdomen

SAFIRE: 2

Sternum Recon*

Sagittal MPR of the Chest

2 mm slice thickness

2 mm slice increment

Kernel: B70s Sharp

Window: Osteo

SAFIRE: 0

Axial Spine Recon*

Without Order

Sagittal MPR entire spine

2 mm slice thickness

2 mm slice increment

Kernel: B70s Sharp

Window: Osteo

SAFIRE: 0

Include the following reconstructions when a Thoracic and Lumbar Spine are ordered

Focused DFOV

Axial Soft Tissue Thoracic Spine

2.5 mm slice thickness

2.5 mm slice increment

DFOV 15 cm

Kernel: I41f medium

Window: Spine

SAFIRE: 2

Axial Bone Thoracic Spine

2.5 mm slice thickness

2.5 mm slice increment

DFOV 15 cm

Kernel: B70s Sharp

Window: Osteo

SAFIRE: 0

Thoracic Spine Sagittal Coronal MPR

1.25 mm slice thickness

1.25 mm slice increment

Kernel: B70s Sharp

Window: Osteo

SAFIRE: 0

Thoracic Spine Coronal MPR

1.25 mm slice thickness

1.25 mm slice increment

Kernel: B70s Sharp

Window: Osteo

SAFIRE: 0

Axial Soft Tissue Lumbar Spine

2.5 mm slice thickness

2.5 mm slice increment

DFOV 15 cm

Kernel: I41f medium

Window: Spine

SAFIRE: 2

Axial Bone Lumbar Spine

2.5 mm slice thickness

2.5 mm slice increment

DFOV 15 cm

Kernel: B70s Sharp

Window: Osteo

SAFIRE: 0

Lumbar Spine Sagittal MPR

1.25 mm slice thickness

1.25 mm slice increment

Kernel: B70s Sharp

Window: Osteo

SAFIRE: 0