

Imaging Protocol – CT Pedi Sinus_Orbits_Facial Bones

PROTOCOL REVIEW 2026 01/df
Approval Imaging Council January 2026

CTDI: 0-4yr: ≤ 15mGy 5-16: ≤ 25mGy

PT Positioning:

- If possible, tilt the patient's head so that a line connecting the lateral canthus of the eye and the EAC is perpendicular to the CT tabletop.

Scan Parameters:

- Sinus - Scan from above frontal sinuses through the bottom of top teeth
- Orbits - Scan from above orbital roof through the maxillary sinuses
- Face - Scan from above orbits through the entire mandible

Setup: LAT scout from no gantry tilt

DFOV: Appropriate to patient body habitus

PACS SERIES

1. SCOUT
2. BONE AX
3. ST AX
4. BONE COR
5. BONE SAG
6. ST COR
7. ST SAG
8. DOSE REPORT

Acquisition Parameters

PT AGE	0-5mo	6-18m	19m-6y	7y-14y	15+
Scan Type	AXIAL	AXIAL	AXIAL	Helical	Helical
Scan FOV	PED HEAD	PED HEAD	SM HEAD	SM HEAD	HEAD
mA *Smart mA/ODM	85	135	220	285	355
kVp *kV assist	100	100	100	100	100
Pitch and Speed (mm/rot)	N/A	N/A	N/A	0.984:1 (140.63 mm/s)	0.984:1 (140.63 mm/s)
Detector Coverage *Smart Coverage	120mm*	140mm*	160mm*	40mm	40mm
Thick	1.25mm	1.25mm	1.25mm	1.25mm	1.25mm
Speed	0.28	0.28	0.28	0.28	0.28

Reconstruction Parameters

RECON 1 (Bone AX)	3000/500
Algorithm	BONE PLUS
ASIR	50%
Recon Type	HELICAL FULL
Slice Thickness	0.625
Increment	0.625
RECON 2 (ST AX)	400/40
Algorithm	STND
ASIR	50%
Recon Type	HELICAL FULL
Slice Thickness	1.25
Increment	1.25
RECON 3 (BONE REF)	3000/500
Algorithm	BONE PLUS
ASIR	50%
Recon Type	HELICAL PLUS
Slice Thickness	0.625
Increment	0.3125
RECON 4 (ST REF)	400/40
Algorithm	STND
ASIR	50%
Recon Type	HELICAL PLUS
Slice Thickness	0.625
Increment	0.3125