

# OtisMed® Imaging Technique: **Quick Reference**

1.5T MRI SAGITTAL



<b>GE FRFSE-XL Sagittal PD</b>	
Time of Scan (TOS)	3-5 minutes
Imaging Options	TRE, Fast, FR, Zip 512      DO NOT USE NPW!!!
Gradient Mode	Whole
Frequency Direction	S-I
Phase Correct	ON
TE (Echo Time)	24
TR (Repetition Time)	3000 (range 2800 – 3200)
ETL (Echo Train Length)	9
NEX	1
Bandwidth	31.25
FOV	16cm (ABSOLUTE VALUE ~ DO NOT CHANGE)
Phase FOV	0 .80 or 0.90 (If large knee use full PFOV)
Slice Thickness	2mm (~50 slices to cover all of the bone, both medial condyle and fibula)
Spacing	0 (Interleaved)
Matrix	256 × 256
Acquisitions	2 – 3
<p><b>3 PLANE LOCATOR:</b></p> <p>Use a 24cm FOV for the 3 Plane locators.</p> <ul style="list-style-type: none"> <li>• 4mm skip 1 mm</li> <li>• Matrix 256 × 224      1 NEX</li> <li>• 9 slices – Axial</li> <li>• 9 slices – Coronal</li> <li>• 5 slices – Sagittal</li> </ul>	

**Philips Achieva 1.5 TSE Sagittal PD**

Time of Scan (TOS)	~ 4 minutes
Imaging Options	Recon 512, 100% Sampling, Asymmetric
Fold Over Direction	AP
TE (Echo Time)	8-13 (Default Echo Spacing between 8-13 seen on info page of scanner)
TR (Repetition Time)	~3000 (range 2800 – 3200)
TSF	8
Watershift	Minimum
Packages	2 – 3 (4 for larger knees)
FOV	160 FH 160 AP 100LR
RFOV	80% to 90% (If large knee use full RFOV)
Slice Thickness	2mm (~50 slices to cover all of bone, both medial condyle and fibula)
Gap	0 or Interleave
Voxel Size	256 × 256 or (FH- - 0.3300 AP - 0.4624)
Recon Matrix	512
Recon Voxel Size	0.3125
NSA	1

If using an older version software try the following:

- Fold Over Suppression – OFF (Make sure phase is A-P direction)

**3 PLANE LOCATOR:**

You may use large FOV for initial transverse locator. Use a 24cm FOV for 3 Plane Locator/Scout scans. Follow with a high resolution 3 Plane Locator/Survey using the following suggested parameters (Approximate Time: 1-2 minutes)

- 4mm skip 1mm
- Matrix 384 × 224      1 NSA
- 9 slices – Axial
- 9 slices – Coronal
- 5 slices – Sagittal

Toshiba 1.5 PD FSE Sagittal	
Pulse Sequence	PD FSE SAGITTAL
Obtain High Res 3 Plane Loc:	<p>You may use large FOV for initial transverse locator. Follow with a high resolution 3 Plane Locator/Survey using a 24cm FOV using the following suggested parameters (Approximate Time: 2 minutes)</p> <ul style="list-style-type: none"> <li>• 4mm skip 1mm</li> <li>• Matrix 256 × 224      1 NEX</li> <li>• 11 slices – Axial</li> <li>• 11 slices – Coronal</li> <li>• 9 slices – Sagittal</li> </ul> <p>DO NOT ANGLE LOCS!!!</p>
Imaging Technique:	FSE+15_sl_t_2mm
Echo Train Spacing (FSE only)	15
Echo Train Length (FSE only)	3
TR (Repetition Time)	~ 3200 (range 3000 - 3600)
TE (Echo Time)	~ 30
Slice Thickness	2
Slice Gap	0
Plane/PE Axis	SG/AP
No Wrap	None
PE Matrix	256 × 256
FOV	16 × 16 (ABSOLUTE VALUE ~ DO NOT CHANGE)
Spatial Preset	None
Imaging Flip Angle	90
Flop Angle	160 (Equal to Tailored RF)
Number of Slices	~ 50
Acquisition Order	Forward
Interleaving	Interleave
Number of Acquisitions (NAQ)	1
Resolution	.62 × .62
Number of Coverages	2
Total Acquisition Time	~ 4 minutes
Filters (Open Recon Tab)	GA00 : Weak
Reconstruction	<p>PE/RO Fine (aka Zip 512)</p> <p>Instructions:</p> <p>To choose PE &amp; RO Fine go to “Sequence Editor”, choose “Open”, choose “Reconstruction” select PE &amp; RO Fine which equals Zip 512.</p>

**Toshiba 1.5 PD FSE Sagittal (Continued)**

**Notes**

To use different images in large window for correct alignment without losing the angle do the following:

- Go to Scan Plan
- Click on “Options”
- Turn “Default Cross” OFF - this holds oblique angle when switching scan planes in main viewer window.

<b>Siemens</b>	
Time of Scan (TOS) min.	~ 4 min
Routine: Slice Group	1
Slices	~ 50
Distance Factor	0
Phase Direction	A to P
Phase Oversampling	50 %
FOV (cm)	160mm (ABSOLUTE VALUE ~ DO NOT CHANGE)
FOV Phase	100
Slice Thickness (mm)	2
TR (Repetition Time)*	2400 (Complete rest of settings and return to Routine Tab and set this value) TR Range 2000-3000 for double concatenation
TE (Echo Time)*	34 (Complete rest of settings and return to Routine Tab and set this value)
Average (NSA)	1
Concatenations	2 or 3
Contrast	Flip Angle = 150 Fat Suppression = None
Resolution (All Filters Off Except)	Elliptical Filter = ON Base Resolution = 256 Phase Resolution = 100 Interpolation box = Checked (Gives reconstruction of 512 matrix)
Geometry	Multi-Slice Mode = Interleaved Series = Interleaved
Sequence: Part 1	Bandwidth = 195 (201) Average Mode = Short-Term
Part 2	Turbo Factor = 7 Pulse Type = Normal Gradient = Fast
Matrix	256 × 256
Misc Tab	ISO
<p><b>Short Bore Scanners</b> – Select Distortion Correction both pre and post processing for Sagittal scan in Knee Coil as well as all Body/Spine Coil Imaging.</p> <p><b>Short Bore Scanners</b> – Select ISO for Sagittal scan in Knee Coil as well as all Body/Spine Coil Imaging.</p>	

**Siemens 1.5 Scanner - 3 Plane Locator**

- Set-up using the first axial locator placing graphics box on slice showing condyles
- Do NOT angle slice planes to knee joint
- Run scan as a higher resolution scan using suggested scan parameters below
- Do not oblique the scan slices to the knee joint. They must be run as true Axial, Sagittal, and Coronal
- Approximate Time: 2 minutes

Slice Group	3 (Axial, Sagittal, Coronal)
Slices	9 (Axial) 5 (Sagittal) 9 (Coronal)
FOV	240mm
TE/TR	As set as default with your scanner
Slice Thickness	4mm 1 skip (25% distance factor)
Averages	2
Resolution	256 × 224

**Notes:**

---

**Joint Replacements**

---

**Trauma, Extremities & Deformities**

---

**Craniomaxillofacial**

---

**Spine**

---

**Biologics**

---

**Surgical Products**

---

**Neuro & ENT**

---

**Interventional Spine**

---

**Navigation**

---

**Endoscopy**

---

**Communications**

---

**Imaging**

---

**Patient Care & Handling Equipment**

---

**EMS Equipment**

325 Corporate Drive  
Mahwah, NJ 07430  
t: 201 831 5000

[www.stryker.com](http://www.stryker.com)

Registration with, and use of, the OtisMed ShapeMatch services by institutions / persons located in the EEA requires that the institutions / persons be registered as being authorized to transfer patient personal and sensitive personal information to the US. By registering with, and using, this service, you represent and warrant that you are in compliance with the data protection laws of your country of residence or registration.

An imaging technologist must always rely on his or her own professional clinical judgment when deciding whether to use a particular CT or MRI scanner on a particular patient. Stryker does not dispense medical advice and expects that imaging technologists are trained in the use of any particular CT or MRI scanner before using it on a particular patient. An imaging technologist must be certified and/or licensed by the appropriate jurisdiction, and/or otherwise qualified to use a CT or MRI scanner, in order to obtain images used to manufacture an OtisMed cutting guide for on a particular patient.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Innovation in Motion, OtisMed, Stryker. All other trademarks are trademarks of their respective owners or holders.

NL10-BR-OM-2539 Rev. 1  
11/11

Copyright © 2011 Stryker  
Printed in USA