

PEDIATRIC ST NECK GE REVOLUTION 256

CTDI: 0-4yr: ≤ 15mGy 5-16: ≤30mGy

PT Prep:

- **IV contrast:**
 - 2cc per kg of 300 mg iodine/non-ionic @ 1-2ml/s, not to exceed 100ml unless otherwise determined by Radiologist.

PT Positioning:

- **Retract shoulders as much as possible & tilt the patient's head so that a line connecting the lateral canthus of the eye & the EAC is perpendicular to the CT tabletop**

Setup: Supine & lateral scout from above the orbital roof to below the aortic arch

DFOV: Appropriate for patient's body habitus.

Scan Parameters:

- **Perform post contrast acquisition from above the orbital roof to below the aortic arch (*with the PT instructed not to swallow*)**
- **Scan delay 60 seconds after start of contrast**

PACS SERIES

1. SCOUT
2. ST AX
3. BONE AX
4. ST CORONAL 1.5X1.5
5. ST SAGITTAL 1.5X1.5
6. DOSE REPORT

PEDIATRIC ST NECK GE REVOLUTION 256

Acquisition Parameters

PT AGE	0-5mo	6-18m	19m-6y	7y-14y	15+
Scan Type	HELICAL	HELICAL	HELICAL	HELICAL	HELICAL
Scan FOV	PED BODY	SM BODY	SM BODY	SM BODY	SM BODY
mA *Smart mA/ODM	10-350*	10-220*	175-450*	10-550*	10-405*
kVp *kV assist	80	100	120	120	120
Focal Spot	S	S	S	S	S
Noise Index	8.2	6.2	9.2	8.2	10.2
Pitch and Speed (mm/rot)	0.984:1 (140.63 mm/s)	0.992:1 (283.48mm/s)	0.992:1 (283.48mm/s)	0.992:1 (158.75 mm/s)	0.992:1 (79.38 mm/s)
Detector Coverage *Smart Coverage	40mm	80	80	80	80
Thick	1.25	1.25	1.25	1.25	1.25
Speed	0.28	0.28	0.28	0.50	1.00

Reconstruction Parameters

RECON 1(Soft Tissue)	
Algorithm	Std
ASIR	50%
Recon Type	Helical Full
Slice Thickness	1.25
Increment	1.25
RECON 2 (Bone)	

PEDIATRIC ST NECK GE REVOLUTION 256

Algorithm	Bone Plus
ASIR	50%
Recon Type	Helical Full
Slice Thickness	1.25
Increment	1.25
RECON 3 (thins for Reformats)	
Algorithm	STND
ASIR	50%
Recon Type	Helical Plus
Slice Thickness	0.625mm
Increment	0.3125mm