

CT DCN CT Cystogram Sets / Foleys

Check which type of foley catheter patient has- You may be able to connect to collection tubing, or directly to foley catheter

- Bard infection control Add-a-foley collection tubing and bag with EZ-Lok
 - Common for older patients with indwelling catheter
 - Connect primary set tubing to EZ-Lok sampling port (luer lock connection)
 - USE STOPCOCK to connect to luer lock on collection tubing
 - Do not need to disconnect bag
 - Clamp collection tubing with hemostats
 - Follow Protocol

- Direct foley catheter connection
 - Temporary catheters placed only for exam, more variety of sizes
 - USE PLASTIC MULTIPURPOSE TUBING ADAPTER to connect primary set tubing to foley (Light or dark blue Christmas tree depending on size)
 - Clamp tubing with hemostats
 - Follow Protocol

Stopcock connection

Tubing adapter

